

**EUR 3494 e**

**EUROPEAN ATOMIC ENERGY COMMUNITY - EURATOM**

**THE DETERMINATION OF CADMIUM IN CADMIUM-TIN ALLOYS**  
**An oscillo-polarographic method**

by

**V. VERDINGH and H. CUYVERS**

**1967**



**Joint Nuclear Research Center  
Geel Establishment - Belgium**

**Central Bureau for Nuclear Measurements - CBNM**



## LEGAL NOTICE

This document was prepared under the sponsorship of the Commission of the European Atomic Energy Community (EURATOM).

Neither the EURATOM Commission, its contractors nor any person acting on their behalf :

Make any warranty or representation, express or implied, with respect to the accuracy, completeness, or usefulness of the information contained in this document, or that the use of any information, apparatus, method, or process disclosed in this document may not infringe privately owned rights; or

Assume any liability with respect to the use of, or for damages resulting from the use of any information, apparatus, method or process disclosed in this document.

This report is on sale at the addresses listed on cover page 4

at the price of FF 2.50

25.—

DM 2.—

Lit. 310

Fl. 1.80

**When ordering, please quote the EUR number and the title, which are indicated on the cover of each report.**

Printed by Smeets  
Brussels, July 1967

This document was reproduced on the basis of the best available copy.



**EUR 3494 e**

THE DETERMINATION OF CADMIUM IN CADMIUM-TIN ALLOYS - An oscillo-polarographic method by V. VERDINGH and H. CUYVERS

European Atomic Energy Community - EURATOM  
Joint Nuclear Research Center - Geel Establishment (Belgium)  
Central Bureau for Nuclear Measurements - CBNM  
Brussels, July 1967 - 12 Pages - 2 Figures - FB 25

A simple polarographic method that allows the determination of Cd in Cd - Sn alloys is described. Sn is removed by treatment of the Cd - Sn solution with a mixture of HBr - Br<sub>2</sub>. Small amounts of Sn reduce the height of the Cd peak. Knowing the Sn peak height a correction can be applied. 1N HCl was chosen as an indifferent electrolyte. The coefficient of variation of the method is less than 1 %.

**EUR 3494 e**

THE DETERMINATION OF CADMIUM IN CADMIUM-TIN ALLOYS - An oscillo-polarographic method by V. VERDINGH and H. CUYVERS

European Atomic Energy Community - EURATOM  
Joint Nuclear Research Center - Geel Establishment (Belgium)  
Central Bureau for Nuclear Measurements - CBNM  
Brussels, July 1967 - 12 Pages - 2 Figures - FB 25

A simple polarographic method that allows the determination of Cd in Cd - Sn alloys is described. Sn is removed by treatment of the Cd - Sn solution with a mixture of HBr - Br<sub>2</sub>. Small amounts of Sn reduce the height of the Cd peak. Knowing the Sn peak height a correction can be applied. 1N HCl was chosen as an indifferent electrolyte. The coefficient of variation of the method is less than 1 %.



**EUR 3494 e**

**EUROPEAN ATOMIC ENERGY COMMUNITY - EURATOM**

**THE DETERMINATION OF CADMIUM IN CADMIUM-TIN ALLOYS**  
**An oscillo-polarographic method**

by

**V. VERDINGH and H. CUYVERS**

**1967**



**Joint Nuclear Research Center  
Geel Establishment - Belgium**

**Central Bureau for Nuclear Measurements - CBNM**

## **SUMMARY**

A simple polarographic method that allows the determination of Cd in Cd - Sn alloys is described. Sn is removed by treatment of the Cd - Sn solution with a mixture of HBr - Br<sub>2</sub>. Small amounts of Sn reduce the height of the Cd peak. Knowing the Sn peak height a correction can be applied. 1N HC was chosen as an indifferent electrolyte. The coefficient of variation of the method is less than 1 %.

# THE DETERMINATION OF CADMIUM IN CADMIUM - TIN ALLOYS

An oscillo-polarographic method<sup>(\*)</sup>

## 1. Introduction

For the analysis of Cd - Sn alloys, used as activation detector for neutron flux measurements, a rapid polarographic method for the determination of Cd was needed.

The Cd content was of the order of 2 - 3 %.

The direct determination of Cd in the presence of large quantities of Sn is impossible because the half-wave potentials are too similar.

( $E_{1/2}$  for Cd in 1M HCl = -0.64 for Sn in 1M HCl - 0.47 S.C.E)  
Sn has to be removed before the determination of Cd.

Sn was volatilized with a HBr - Br<sub>2</sub> mixture recommended by [2].

The Cd determination itself can be done in a variety of electrolytes [1], [4].

## 2. Experimental

Apparatus: A Southern A 1660 Davis Differential Cathode Ray polarograph was used. The potentials were measured versus a mercury pool anode. Some modifications of the apparatus are described in a previous publication [3].

### 2.1. Reagents

Cd : hyper pure : 99.999 % U.C.B. nr. 5417

Sn : hyper pure : 99.999 % U.C.B.

HCl Merck pro analysi

Merck titrisol

HBr Merck pro analysi

Br<sub>2</sub> Merck pro analysi

HNO<sub>3</sub> Merck supra pure

H<sub>2</sub>O doubly distilled

(\*) Manuscript received on May 23, 1967.

## 2.2. Preparation of solutions

### 2.2.1. Cadmium standard solution

About 500 mg Cd were dissolved in 100 ml conc. HCl + 1 ml conc.  $\text{HNO}_3$ . The solution was quantitatively transferred to a weighed 500 ml flask, diluted and weighed.

From this solution, reference solutions were prepared by weighing. Corrections for buoyancy were made for all weighing operations.

### 2.2.2. Sn solution

About 5.5 g Sn were dissolved in 20 ml conc.  $\text{HNO}_3$  + 100 ml conc. HCl; after dissolution 40 ml conc. HCl were added and the volume was made up to 1000 ml in a volumetric flask.

### 2.2.3. HBr - Br mixture

50 ml HBr 40 % and 10 ml  $\text{Br}_2$  were mixed in a 100 ml flask.

### 2.2.4. Cadmium solutions

Weighed portions of about 10 - 12 g Cd-stock solution were transferred to 250 ml beakers. For the Cd - Sn solutions, 30 ml Sn-stock solution was added.

After treatment with a HBr -  $\text{Br}_2$  mixture, the residue was dissolved in 1 N HCl and quantitatively transferred to 100 ml flasks, diluted with 1 N HCl and weighed.

### 2.2.5. Sample solutions

About 180 mg of the Cd - Sn alloy were dissolved in 5 ml conc. HCl and 1 ml conc.  $\text{HNO}_3$ . After treatment with the HBr -  $\text{Br}_2$  mixture the residue was dissolved



in 1 N HCl, quantitatively transferred to a 100 ml flask, diluted with 1 N HCl and weighed.

### 2.3. Removing of Sn

Cd can be determined in the presence of 2 % of tin. In the presence of larger quantities of Sn the Sn has to be removed until a concentration smaller than 10 % is reached.

A correction can be applied for the influence of the remaining Sn on the Cd peak (see 3.2.).

The Cd - Sn solutions were treated three times with 5 ml HBr - Br<sub>2</sub> mixture and carefully evaporated to dryness on a sand bath at a temperature of about 150°C.

The Sn is removed nearly quantitatively.

The influence of the remaining Sn on the Cd peak height after three evaporations is less than 1 % on the Cd content. A correction for this effect can be made. (see paragraph 3.2.).

### 2.4. Method

5 ml of the sample or calibration solution were brought into a cell and deaerated for 10 minutes with high purity nitrogen. "Slope compensation" and "condensor compensation" were used.

The start potential was about - 0.380 V. The peak potential for Cd lies at about - 0.670 Volt versus Hg-pool anode. In order to improve the reproducibility of the method, at least four or five cell-fillings were made from each solution. Each cell-filling was measured 11 times.

The mean of these measurements was used for further calculations.

### 3. Results

#### 3.1. Calibration factor for Cd (C.F.)

The calibration factor for Cd is derived from measurements of pure Cd-reference solutions without treatment with the HBr - Br<sub>2</sub> mixture. It is the number of scale divisions at maximum sensitivity divided by the concentration of Cd given in µg/g solution.

$$\text{C.F.} = 1680 \text{ S.D. (scale divisions)} / \mu\text{g Cd/g solution}$$

This value is used as a reference throughout the whole study. All peak heights were normalised to this value, even when changing capillary characteristics. The calibration factor of Sn is derived from measured Sn peak heights in presence of Cd (see 3.2.).

$$\text{C.F.} \approx 380 \text{ S.D.} / \mu\text{g Sn/g solution}$$

At concentrations of less than 10 % Sn (with respect to the Sn - Cd content) the reproducibility on the Sn peak height reading is about 20 %. Because the correction factor for influence of Sn on the Cd peak height is less than 1.010, a reproducibility of 20 % on this factor lies within the reproducibility of the method.

#### 3.2. Influence of Sn on the Cd peak height

Solutions with known concentration ratio of Cd and Sn were measured. The decrease of the Cd peak height is proportional to the height of the Sn peak.



Calc.concen- tration in ug/g sol.		% Sn in Cd+Sn	Peak height (S.D.)		Found conc. ugCd/g sol.	Deviation from the calc.value (ug Cd/g sol.)
Cd	Sn		Cd	Sn		
80.43	0	0	135,100	0	80.43	- 0.00
80.43	9.10	10.20	134,300	3,500	79.94	- 0.49
80.43	22.75	22.00	130,400	8,700	77.62	- 2.81
80.43	45.50	36.10	125,800	17,500	74.88	- 5.55

A correction factor on the Cd peak height, proportional to the Sn peak height, can be applied.

Sn peak height S.D.	Correction factor for the Cd peak height
0	1.000
3,500	1.006
8,700	1.036
17,500	1.074

These values are illustrated in fig. 1. By extrapolation of the curve for lower Sn concentrations (after treatment with the HBr - B<sub>2</sub> mixture) one obtains another curve, fig. 2., from which the correction factor can be read. Working under normal conditions the correction factor to be applied is smaller than 1.010.

### 3.3. Analysis of Cd - Sn solutions after treatment with the HBr - Br<sub>2</sub> mixture

In all solutions the remaining Sn peak height was about 1500 S.D. ( $\pm 20\%$ ). The correction factor for 1500 S.D. Sn is 1.0062.

Peak height Cd(SD) measured	Corrected peak height Cd(SD)	Calc.conc. Cd in $\mu\text{g/g}$ solution	Found conc. Cd in $\mu\text{g/g}$ solution	Deviation from the calc.value $\mu\text{gCd/gsol.}$
81,400	81,950	48.67	48.76	+ 0.10
81,600	82,150	48.97	48.88	- 0.09
84,350	84,900	50.79	50.53	- 0.26
72,550	73,000	43.87	43.51	- 0.36
82,300	82,850	49.21	49.31	+ 0.10
85,050	85,600	50.59	50.94	+ 0.35
79,500	80,000	47.32	47.68	+ 0.36
40,700	91,300	53.73	54.35	+ 0.62
84,300	84,850	50.67	50.18	- 0.49
84,200	84,200	50.65	50.25	- 0.40
72,500	72,500	43.19	43.19	0.00

The last two solutions were pure Cd solutions without Sn, but also treated with the HBr - Br mixture.

The coefficient of variation for this series of measurements  $S/\bar{x} = 0.61\%$ .



### 3.4. Analysis of Cd - Sn alloys

From each sample three solutions were made and each solution was measured 3 times. In all solutions the remaining Sn peak height was of the order of 1000 S.D. The correction factor for 1000 S.D. Sn is 1.0042 (see 3.2.).

Sample	Cd peak height S.D.meas.	Cd peak height S.D.corr.	Weight sample g	Weight sol. g	% Cd in Cd-Sn
Tp 326  R	87,450	87,850	0.19835	102.0410	2.731
	87,900	88,300			2.744
	88,000	88,400			2.747
	85,050	85,450	0.19215	101.6329	2.730
	85,000	85,400			2.729
	84,700	85,100			2.719
	82,500	82,850	0.18830	102.5136	2.726
	82,650	83,000			2.731
	82,600	82,950			2.730
$\bar{x} = 2.732 \qquad S/\bar{x} = 0.32 \%$					
Tp 326  M	79,050	79,400	0.17805	101.4835	2.737
	79,200	79,550			2.742
	76,900	79,250			2.732
	85,300	85,700	0.19322	101.5378	2.721
	85,300	85,700			2.721
	85,500	85,900			2.727
	82,700	83,050	0.18540	101.4326	2.739
	82,500	82,850			2.732
	83,000	83,350			2.748
$\bar{x} = 2.733 \qquad S/\bar{x} = 0.34 \%$					

### Acknowledgements

The authors wish to express their gratitude to Dr. K.F. LAUER for his interest and support.

#### 4. References

- [1] I.M. Kolthoff and Lingane: Polarography New York Intersci. Publ.
- [2] A.S.T.M. Methods of Chemical Analysis of Metals 1946, p. 239
- [3] V. Verdingh, Z. Anal. Chem. Bd. 220 Band 3, Heft 1966.
- [4] V. Verdingh - M. Van Bouwel - unpublished.






fig.1




fig.2






## NOTICE TO THE READER

All Euratom reports are announced, as and when they are issued, in the monthly periodical **EURATOM INFORMATION**, edited by the Centre for Information and Documentation (CID). For subscription (1 year : US\$ 15, £ 5.7) or free specimen copies please write to :

**Handelsblatt GmbH**  
**"Euratom Information"**  
**Postfach 1102**  
**D-4 Düsseldorf (Germany)**

OR

**Office central de vente des publications**  
**des Communautés européennes**  
**2, Place de Metz**  
**Luxembourg**



To disseminate knowledge is to disseminate prosperity — I mean general prosperity and not individual riches — and with prosperity disappears the greater part of the evil which is our heritage from darker times.

Alfred Nobel



## SALES OFFICES

All Euratom reports are on sale at the offices listed below, at the prices given on the back of the front cover (when ordering, specify clearly the EUR number and the title of the report, which are shown on the front cover).

### OFFICE CENTRAL DE VENTE DES PUBLICATIONS DES COMMUNAUTES EUROPEENNES

2, place de Metz, Luxembourg (Compte chèque postal N° 191-90)

#### BELGIQUE — BELGIË

MONITEUR BELGE  
40-42, rue de Louvain - Bruxelles  
BELGISCH STAATSBAD  
Leuvenseweg 40-42, - Brussel

#### LUXEMBOURG

OFFICE CENTRAL DE VENTE  
DES PUBLICATIONS DES  
COMMUNAUTES EUROPEENNES  
9, rue Goethe - Luxembourg

#### DEUTSCHLAND

BUNDESANZEIGER  
Postfach - Köln 1

#### NEDERLAND

STAATSDRUKKERIJ  
Christoffel Plantijnstraat - Den Haag

#### FRANCE

SERVICE DE VENTE EN FRANCE  
DES PUBLICATIONS DES  
COMMUNAUTES EUROPEENNES  
26, rue Desaix - Paris 15<sup>e</sup>

#### UNITED KINGDOM

H. M. STATIONERY OFFICE  
P. O. Box 569 - London S.E.1

#### ITALIA

LIBRERIA DELLO STATO  
Piazza G. Verdi, 10 - Roma

EURATOM — C.I.D.  
51-53, rue Belliard  
Bruxelles (Belgique)

CDNA03494ENC