

EUROPEAN ATOMIC ENERGY COMMUNITY - EURATOM

# THE DETERMINATION OF FREE ACIDITY IN PLUTONIUM CONTAINING SOLUTIONS AND THE SEMI-QUANTITATIVE ACIDIMETRIC ESTIMATION OF PLUTONIUM

by

化加引出的计算管用

K. BUIJS, B. CHAVANE DE DALMASSY and M.J. MAURICE

1968



Joint Nuclear Research Center Karlsruche Establishment — Germany

European Institute for Transuranium Elements

## LEGAL NOTICE

This document was prepared under the sponsorship of the Commission of the European Communities.

Neither the Commission of the European Communities, its contractors nor any person acting on their behalf:

Make any warranty or representation, express or implied, with respect to the accuracy, completeness, or usefulness of the information contained in this document, or that the use of any information, apparatus, method, or process disclosed in this document may not infringe privately owned rights; or

Assume any liability with respect to the use of, or for damages resulting from the use of any information, apparatus, method or process disclosed in this document.

This report is on sale at the addresses listed on cover page 4

at the price of FF 2.50	FB 25.—	DM 2.—	Lit. 310	Fl. 1.80
	The set of the set of the set of the set of the	A SHARE SHEET ALL ALL ALL ALL ALL ALL ALL ALL ALL AL	A A A A A A A A A A A A A A A A A A A	I CIMER I LINE MARKED I I

When ordering, please quote the EUR number and the title, which are indicated on the cover of each report.

Printed by Van Muysewinkel Brussels, June 1968.

This document was reproduced on the basis of the best available copy.

#### EUR 3934 e

THE DETERMINATION OF FREE ACIDITY IN PLUTONIUM CON-TAINING SOLUTIONS AND THE SEMI-QUANTITATIVE ACIDIMETRIC ESTIMATION OF PLUTONIUM by K. BUIJS, B. CHAVANE DE DALMASSY and M.J. MAURICE

European Atomic Energy Community — EURATOM Joint Nuclear Research Center — Karlsruhe Establishment (Germany) European Institute for Transuranium Elements Brussels, June 1968 — 10 Pages — FB 25

A fast method for the analysis of free acid in plutonium solutions has been developed. An ammonium fluoride solution provides the medium in which hydrolysis of plutonium is prevented.

#### EUR 3934 e

THE DETERMINATION OF FREE ACIDITY IN PLUTONIUM CON-TAINING SOLUTIONS AND THE SEMI-QUANTITATIVE ACIDIMETRIC ESTIMATION OF PLUTONIUM by K. BUIJS, B. CHAVANE DE DALMASSY and M.J. MAURICE

European Atomic Energy Community — EURATOM Joint Nuclear Research Center — Karlsruhe Establishment (Germany) European Institute for Transuranium Elements Brussels, June 1968 — 10 Pages — FB 25

A fast method for the analysis of free acid in plutonium solutions has been developed. An ammonium fluoride solution provides the medium in which hydrolysis of plutonium is prevented.

\_\_\_\_\_

····

# EUR 3934 e

### EUROPEAN ATOMIC ENERGY COMMUNITY - EURATOM

# THE DETERMINATION OF FREE ACIDITY IN PLUTONIUM CONTAINING SOLUTIONS AND THE SEMI-QUANTITATIVE ACIDIMETRIC ESTIMATION OF PLUTONIUM

by

K. BUIJS, B. CHAVANE DE DALMASSY and M.J. MAURICE

1968



Joint Nuclear Research Center Karlsruche Establishment — Germany

European Institute for Transuranium Elements

### Summary

A fast method for the analysis of free acid in plutonium solutions has been developed. An ammonium fluoride solution provides the medium in which hydrolysis of plutonium is prevented.

### **KEYWORDS**

PLUTONIUM COMPOUNDS SOLUTIONS ACIDITY DETERMINATION AMMONIUM COMPOUNDS FLUORIDES STABILITY HYDROLYSIS

#### · ·

## CONTENTS

Introduction	5	
Experimental part	6	
Discussion	. 8	/ 1
Conclusions ,	10	
References	10	

.

• . .

The determination of free acidity in plutonium containing solutions

and the semi-quantitative acidimetric estimation of plutonium (+)

#### INTRODUCTION

The analysis of free acid in plutonium solutions, as in all heavy metal solutions, is hampered by the hydrolysis of the plutonium ions, which occurs at a pH of about 1 - 2. If no precautions are taken, the acidity values found are in-evitably high.

Various procedures have been developed to eliminate excess alkali consumption due to hydrolysis. Methods most frequently used are elimination of plutonium from the solution by precipitation prior to analysis or complexation of the metal. Potassium iodate is often used for precipitation (1). It is known that this method introduces a negative bias. Moreover, it involves a time-consuming separation step and causes difficulties with the concentrated iodate solutions during waste treatment. For complexation, a solution of sodium citrate (2) has been used. This procedure, however, introduces a positive bias due to the fact that the plutonium citrate complex is not stable over the entire pH range of the titration. The resulting waste solution still requires special attention.

In view of this situation we decided to study a new method for the free acid analysis, which has to be accurate, fast, without separation step, and leads to acceptable waste solutions.

(+) Manuscript received on March 24, 1968.

#### Experimental part

The complexing or precipitating agent should be a neutral salt or a salt of a relatively strong acid and a weak base, to avoid any positive bias. The plutonium complex or precipitate sould be stable over the whole pH range of the titration. A series of experiments were conducted with ammonium fluoride which fulfills these requirements.

On solutions of plutonium nitrate in nitric acid three determinations were performed :

- a) An aliquot of the sample solution (1 ml for acidities higher than 4 N or 2 ml for acidities lower than 4 N) was added to 10 ml of 1.5 M ammonium fluoride solution. The solution was titrated potentiometrically with 1 N sodium hydroxide solution to the original pH of the ammonium fluoride solution with disregard of the plutonium fluoride precipitate which was formed. A glass electrode and a saturated calomel electrode were used. In this way, the free acidity, x, (see table 1), was determined.
- b) For each sample solution a determination of total acidity, y (table 1), was also made in the absence of fluoride by potentiometric titration of a suitable aliquot with 1 N sodium hydroxide solution.
- c) In each sample solution the plutonium content, z, was determined according to Corpel and Regnaud (3).

When all results are expressed in mequivs, the following equation holds :

$$\mathbf{x} = \mathbf{y} - \mathbf{z} \tag{1}$$

The results are given in table 1.

- 6 -

<u>Table 1</u>

free acid		total acid	plutonium		free acid	difference		
sample	x	У	<b>Z</b>	Z '	y - 2	$\Delta = x - (y - z)$		
1	1.74 N	2.00 N	0.26 N	15.3g/L	1.74 N	ON		
2	3.34	3.65	0.26	15.3	3.39	-0.05		
3	1.36	1.73	0.39	23.2	1.34	+0.02		
4	1.22	1.83	0.52	31.0	1.31	-0.09		
5	3.67	4.15	0.54	32.1	3.61	+0.06		
6	1.50	2.10	0.58	34.8	1.52	+0.02		
7	1.85	2.65	0.65	38.8	2.00	-0.15		
8	1.34	1.97	0.66	39•4	1.31	+0.03		
9	1.14	1.80	0.72	42.7	1.08	+0.06		
10	2.30	2.93	0.72	42.7	2.21	+0.09		
11	1.74	2.37	0.73	43.5	1.64	+0.10		
12	0.92	1.77	0.86	51.2	0.91	+0.01		
average, $\overline{\Delta}$ = +0.012 N								
standard deviation of $\Delta$ , s $\Delta = 0.07$ N								

.

.

#### Discussion

It cannot be detected that the  $\Delta$ -values of table 1 contain outliers and, therefore, all values can be used for further calculations. The average difference between the results of the two types of measurements,  $\overline{\Delta}$ , does not differ significantly from zero on the 0.05 level, which follows from application of the t-test. This means that the method has no detectable bias. A correlation between the  $\Delta_i$ -values and the plutonium concentrations,  $z_i$ , cannot be detected with a modified sign test.

From 
$$\Delta = x - y + z$$
 (II)  
it follows  $s_{\Delta}^2 = s_x^2 + s_y^2 + s_z^2$  (III)  
where:  $s_x = standard$  deviation of a single free  
acidity determination  
 $s_y = standard$  deviation of a single total  
acidity determination  
 $s_z = standard$  deviation of a single pluto-  
nium determination.

Since  $s_z^2$  is negligible with respect to  $s_x^2$  and  $s_y^2$ , and moreover,  $s_x^2$  may be put equal to  $s_y^2$ , it is seen that

$$s_{\Delta}^{2} = 2s_{x}^{2}$$
or  $s_{x} = 1/2 s_{\Delta}\sqrt{2}$ 
(IV)

Substitution of s = 0.07 in (IV) gives :  

$$s_x = s_y = 0.05$$
 (V)

As results are unbiased, it can be said that

 $z = y - x \tag{VI}$ 

This means that in a solution of pure plutonium nitrate in nitric acid the plutonium concentration can be estimated very rapidly by carrying out an acidity determination in the absence and one in the presence of ammonium fluoride.

From (VI) it follows, that

$$s_z^2 = s_y^2 + s_x^2 = 2 s_x^2$$
 (VII)

Substitution of  $s_x = 0.05$  in (VII) gives

$$s_z = 0.07$$
 equivs/L (VIII)  
or  $s_z' = 4.2$  g/L

From eq.(VIII) is follows that this plutonium determination is essentially a semi-quantitative estimation.

#### CONCLUSIONS

- The method gives unbiased results for solutions containing
   15 to 50 g of plutonium/liter and having acidities between
   0.9 and 4 N.
- 2. The standard deviation of a single determination of the free acid concentration is 0.05 N.
- 3. Since results are unbiased (conclusion 1), plutonium can be estimated in a solution of its nitraté in nitric acid as the difference between total acidity and free acidity.
- 4. The standard deviation of a single plutonium estimation, mentioned under 3, is 4.2 g/L, which means that this estimation is a semi-quantitative one.

#### References

- 1) W. BÄHR and D. THIELE, Report K.F.K. 499 (1966).
- 2) J.L. PFLUG and F.J. MINER, Anal. Chim. Acta, 23, 362 (1960).
- 3) J. CORPEL and F. REGNAUD, Anal. Chim. Acta, 35, 508 (1966).

## NOTICE TO THE READER

ANH HITS

All Euratom reports are announced, as and when they are issued, in the monthly periodical **EURATOM INFORMATION**, edited by the Centre for Information and Documentation (CID). For subscription (1 year : US\$ 15,  $\pm$  6.5) or free specimen copies please write to :

Handelsblatt GmbH "Euratom Information" Postfach 1102 D-4 Düsseldorf (Germany)

or

Office central de vente des publications des Communautés européennes 2, Place de Metz Luxembourg

To disseminate knowledge is to disseminate prosperity — I mean general prosperity and not individual riches — and with prosperity disappears the greater part of the evil which is our heritage from darker times.

Alfred Nobel

# SALES OFFICES

All Euratom reports are on sale at the offices listed below, at the prices given on the back of the front cover (when ordering, specify clearly the EUR number and the title of the report, which are shown on the front cover).

OFFICE CENTRAL DE VENTE DES PUBLICATIONS DES COMMUNAUTES EUROPEENNES

2, place de Metz, Luxembourg (Compte chèque postal Nº 191-90)

BELGIQUE — BELGIË MONITEUR BELGE 40-42, rue de Louvain - Bruxelles BELGISCH STAATSBLAD Leuvenseweg 40-42 - Brussel

DEUTSCHLAND BUNDESANZEIGER Postfach - Köln 1

FRANCE

SERVICE DE VENTE EN FRANCE DES PUBLICATIONS DES COMMUNAUTES EUROPEENNES 26, rue Desaix - Paris 15°

ITALIA LIBRERIA DELLO STATO Piazza G. Verdi, 10 - Roma LUXEMBOURG OFFICE CENTRAL DE VENTE DES PUBLICATIONS DES COMMUNAUTES EUROPEENNES 9, rue Goethe - Luxembourg

NEDERLAND STAATSDRUKKERIJ Christoffel Plantijnstraat - Den Haag

UNITED KINGDOM H. M. STATIONERY OFFICE P. O. Box 569 - London S.E.1



Į

CDNA03934ENC