

COMMISSION OF THE EUROPEAN COMMUNITIES

A DIGITAL SYSTEM FOR SPEED AND PHASE CONTROL OF PHASED NEUTRON CHOPPER ARRAYS

by

E. BETTENDROFFER and J. EDER

1974



Joint Nuclear Research Centre Ispra Establishment - Italy

Electronics Division

LEGAL NOTICE

This document was prepared under the sponsorship of the Commission of the European Communities.

Neither the Commission of the European Communities, its contractors nor any person acting on their behalf:

make any warranty or representation, express or implied, with respect to the accuracy, completeness, or usefulness of the information contained in this document, or that the use of any information, apparatus, method or process disclosed in this document may not infringe privately owned rights; or

assume any liability with respect to the use of, or for damages resulting from the use of any information, apparatus, method or process disclosed in this document.

This report is on sale at the addresses listed on cover page 4

at the price of B.Fr. 40.-

Commission of the European Communities D.G. XIII - C.I.D. 29, rue Aldringen Luxembourg

May 1974

This document was reproduced on the basis of the best available copy

EUR 5112 e

1

!

1

1

1

1

A DIGITAL SYSTEM FOR SPEED AND PHASE CONTROL OF PHASED NEUTRON CHOPPER ARRAYS by E. BETTENDROFFER and J. EDER

Commission of the European Communities Joint Nuclear Research Centre - Ispra Establishment (Italy) Electronics Division Luxembourg, May 1974 - 28 Pages - 11 Figures - B.Fr. 40,—

Phased chopper arrays for slow neutron experiments require electrical drive systems with high stability in speed and phasing of the different rotors. In (1) a modular drive system is described which uses digital speed and phase control. This report gives the design of the control modules for speed and phase settings.

All timing signals are derived from a central crystal clock which guarantees the required high longterm stability of speeds and phases. The settings may be adjusted manually or automatically by a computer through an appropriate interface (2).

EUR 5112 e

A DIGITAL SYSTEM FOR SPEED AND PHASE CONTROL OF PHASED NEUTRON CHOPPER ARRAYS by E. BETTENDROFFER and J. EDER

Commission of the European Communities Joint Nuclear Research Centre - Ispra Establishment (Italy) Electronics Division Luxembourg, May 1974 - 28 Pages - 11 Figures - B.Fr. 40,—

Phased chopper arrays for slow neutron experiments require electrical drive systems with high stability in speed and phasing of the different rotors. In (1) a modular drive system is described which uses digital speed and phase control. This report gives the design of the control modules for speed and phase settings.

All timing signals are derived from a central crystal clock which guarantees the required high longterm stability of speeds and phases. The settings may be adjusted manually or automatically by a computer through an appropriate interface (2).

EUR 5112 e

A DIGITAL SYSTEM FOR SPEED AND PHASE CONTROL OF PHASED NEUTRON CHOPPER ARRAYS by E. BETTENDROFFER and J. EDER

Commission of the European Communities Joint Nuclear Research Centre - Ispra Establishment (Italy) Electronics Division Luxembourg, May 1974 - 28 Pages - 11 Figures - B.Fr. 40,—

Phased chopper arrays for slow neutron experiments require electrical drive systems with high stability in speed and phasing of the different rotors. In (1) a modular drive system is described which uses digital speed and phase control. This report gives the design of the control modules for speed and phase settings.

All timing signals are derived from a central crystal clock which guarantees the required high longterm stability of speeds and phases. The settings may be adjusted manually or automatically by a computer through an appropriate interface (2).

1 . .

EUR 5112 e

COMMISSION OF THE EUROPEAN COMMUNITIES

A DIGITAL SYSTEM FOR SPEED AND PHASE CONTROL OF PHASED NEUTRON CHOPPER ARRAYS

by

E. BETTENDROFFER and J. EDER

.

1974



Joint Nuclear Research Centre Ispra Establishment - Italy

Electronics Division

ABSTRACT

Phased chopper arrays for slow neutron experiments require electrical drive systems with high stability in speed and phasing of the different rotors. In (1) a modular drive system is described which uses digital speed and phase control. This report gives the design of the control modules for speed and phase settings.

All timing signals are derived from a central crystal clock which guarantees the required high longterm stability of speeds and phases. The settings may be adjusted manually or automatically by a computer through an appropriate interface (2).

CONTENTS

•

1)	PRINCIPLE OF OPERATION	5					
2)) MASTER GENERATOR						
	2.1 Manual-mode	8					
	2.2 Automatic-mode	10					
	2.3 Specifications	11					
3)	DELAY GENERATOR	13					
	3.1 Manual-mode	13					
	3.2 Automatic-mode	14					
	3.3 Specifications	15					
4)	BIBLIOGRAPHY	16					

1) PRINCIPLE OF OPERATION

Fig. 2 shows the control system for an array of 4 phased choppers. Each chopper has its own power drive module which contains also the control logic for motor start-up and positional feedback which forces the chopper to remain phaselocked with the timing signals F and 6F.

A master generator provides the timing signals for the first chopper and determines the speed of the chopper array. It delivers also a pilot signal which serves as a general time base for all chopper control units. Phase angles between different rotors are always referred to this pilot. Delay generators produce by digitally delaying the pilot pulse frame the timing signals for the delayed choppers with the necessary frequency and phase relation.

With this arrangement the speed of all choppers can be changed or adjusted simultaneously, while the phase relations can be adjusted individually.

The timing signals F and 6F are pulse trains with pulse frequencies of one pulse resp. six pulses per chopper revolution. The period time T of the 6F signal is used for speed setting.

- 5 -

The relation between chopper speed and period T is given by :

$$\omega = \frac{1}{6.T} \cdot 60 = 10^7 \cdot \frac{1}{T} (\mu s) rpm$$
 (1)

T can be chosen between

in steps of 1 jus. This corresponds to a resolution of 0.5 % at 25000 rpm resp. 0.02 % at 1000 rpm. The speed stability however depends only on the crystal clock and is in the order of 5.10^{-4} .

In a similar manner instead of phase angles the delay time D between the aperture times of a chopper referred to the master chopper is used as set value for the delay generators due to

$$\varphi = 360. \ \underline{\omega} \cdot D = 60.10^{-6} \cdot \omega \cdot D_{(yus)} \checkmark^{\circ} (2)$$

Since the delay can be adjusted in steps of 1 us too, the angular resolution is

$$0.06^{\circ} = 3.6'$$
 at 1000 rpm
and
 1.5° at 25000 rpm

Manual or automatic (via computer interface) input of the set values for speed and delay can be chosen.

Since the chopper-drives are phaselocked to the timing signals F and 6F, any step change of the signal frequency must be avoided if we tant to keep all choppers in synchronism. Therefore a special circuit is built in, which smoothes the steps provoked by an arbitrary change of a set value in a way, the phase lock system can follow without losing synchronism. For safety reasons the maximum speed must be limited to ≤ 25000 rpm. This corresponds to a period set value of

≱ 400 /us. If by error a period value below this limit is selected (manually or automatically), the frequency of the timing signals is automatically limited.

and the second second

2) MASTER GENERATOR (Fig. 4)

2.1. Manual-mode

When supply voltage is applied, the one-shots I and II are triggered by "RC" (fig. 5). One-shot I causes the preselected number (3000 for example) to be stored in the latches (fig. 6). One-shot II produces a parallel load of the up'down counters, which means that the value 3000 is fed to the counter outputs through the data input line. Since the data inputs of the 0.1 /us decade (fig. 7) are at logical 0, this counter is loaded with zero. Thus, in reality, the loaded value is not 3000 but 30000. The A and B inputs of the upper comparators both "see" the number 30000, their outputs A > B and A < B are at logical zero, and the output pulses of the slew generator (fig. 5)are inhibited.

The 400 /us interlock is controlled by the 400 /us coincidence : When the contents of the up/down counters become smaller than 400 /us, the outputs of the interlock (representing decades 1, 2 and 3) remain at 400 (fig. 6). In our example (P = 3000 /us) the value 000 is fed through the interlock and the A inputs of the lower comparators are 30000. When the synchronous decade counter (fig. 7), supplied by the 10 Mhz crystal oscillator, reaches 30000, the A > B outputs and consequently the output of the Or-gate (fig. 8) become zero. With this edge the decade counter is reset (and recommences to count) and one-shot III is triggered. This one-shot gives the 6F output signal (after level translation) and the clock pulses for the divide by six counter, which will deliver an output pulse at each sixth 6F pulse. This pulse (after level translation) is used as 1 F output signal and also resets the divide by six counter to zero.

Assume the value 4000 is now read into the latches by pressing the load pushbutton. Now, the A > B output of the upper comparators becomes "1" and the slew generator pulses are fed to the up input of the up/down counters until they have reached 40000.

After the first slew generator pulse, the up/down counter outputs and also the A inputs of the lower comparators are 30001. Thus, the decade counter will count up to 30001 before the lower comparators deliver an output signal. In this way the period of the 6F output signal follows the up/down counter step by step, increasing by steps of 0.1 /us at a rate depending on the slew generator frequency. The decreasing mechanism is similar.

The slew generator operates as follows : The oscillator is an unijunction transistor with recurrence frequency depending on the RC value. The fourth decade and the three last bits of the third decade of the up/down counters are used to produce a coincidence at 4 different period values. Each coincidence commutates to ground a capacitor at predeterminated period values, thus increasing or decreasing the total capacitor value and reciprocal the slew rate.

Ī	Perio	od	val	ue		Slew	rate
9999	us	≥	Р ≥	4000	us	30	ms
4000	11	>	P ≥	2000	11	70	ms
2000	11	>	Ρ≥	1000	17	200	ms
1000	"	>	₽≥	600	17	300	ms
600	**	≥	P ≥	400	19	600	ms

2.2. Automatic-mode

The period value to be stored into the latches is supplied by a computer, which also delivers the load signal. After the value has been inscribed, the generator delivers an acknowledge pulse to the computer. See also [2].

2.3. Specifications

The Master Generator is housed in a two-units ESONE plug-in (see fig. 1). A second plug-in contains the common power supply for one Master Generator and up to three Delay Generators. <u>Note</u> : The generators and their power supply should be operated

only in the supplied chassis.

Since the standard ESONE wiring has been modified at the present device, it would also be dangerous to operate it in the chassis instruments which are designed in the ESONE standard.

Don't exchange the position of the power supply with a generator and vice versa. However, the positions of the generators may be chosen arbitrarily.

6 F output signal : period P1 variable between 400 and 9999 us, by steps of 0.1 us.

1 F output signal : period P2 = 6P1, synchronous with 6F (fig. 3). Variation : 2400 to 59994 /us by steps of 0.6 /us.

Common characteristics of both output signals :

- pulse width : 2 μ s \pm 10 %
- quiescent level : + 15 V
- pulses negative from + 15 V to ground
- capability of driving up to 10 MHTL circuits (I sink = 30 mA) from each output connector, including the test outputs.

- slew rate : \simeq 30 ms to \simeq 600 ms in 5 ranges. Automatic range commutation at predetermined values of P1. The rate value can be adjusted by a factor 1 : 3 by front panel potentiometer. the second se

- Readout : instantaneous value of P1 (four decades) displayed by seven segment indicators.
- 10 Mhz, produced by a crystal oscillator Clock frequency :
- Operating modes : Manual : presetting and loading of the desired period value are executed by the operator.
 - Automatic : the functions are computer controlled.

Power supply : + 5 V, 10 A (N.V. Delta Elektronika Type M5-10) + 15 V, 150 mA.

3) DELAY GENERATOR (Fig. 9)

3.1. Manual-mode

For the operation of the decade command card (fig. 6) and the command card (fig. 5) refer to 2.1. of the Master Generator.

The 1F pulse issued from the Master Generator starts the delay counter (fig. 10), supplied by the 1 us clock pulses from Master Generator. When the counter reaches the presetted delay value, the output pulse of the comparator simultaneously : 1) resets the delay counter and inhibits the clock pulses 2) resets the D > 6P interlock

3) triggers the 1F' output one-shot.

The D > 6P interlock operates as follows : When the second 1F pulse comes in before the comparator delivers an output pulse (that means : when the presetted delay value > 1F period), this 1F pulse is fed through the interlock and takes over the 3 functions described above ; during this time, the comparator pulse is inhibited. When the delay value becomes \leq 6P1, the comparator pulse passes, while the 1F pulse is inhibited.

When D = 0, the output of the comparator is permanently logical "1" and the 1F pulse is used as output signal.

1F' resets a divide by six counter, the outputs of which control the select inputs of a multiplexer (fig. 11).

The same signal, slightly delayed is applied to the corresponding data input of the multiplexer and appears at its output and triggers the one-shot. The latter delivers the output signal 6F', the trailing edge of which is used to advance the divide by six counter by one position, enabling one of the five paralleled data inputs of the multiplexer.

The multiplexer output signal also resets a decade counter, which is supplied with the clock pulse issued from PC20224 (Master Generator). The contents of this decade counter are compared with the value of P1 (Master Generator); when this value is reached, the comparator output signal is applied to one of the 5 paralleled inputs of the multiplexer. This process is repeated five times. The fifth output pulse sets the divide by six counter to a position corresponding to zero, so enabling the 1F' signal. In this way, each 1F' pulse is used as reference to the 6F' pulses, thus adjusting 6F' to the instantaneous delay value.

3.2. Automatic-mode

The delay value and the "load" pulse are supplied by the computer (same procedure as with Master Generator, see 2.2.).

- 14 -

3.3. <u>Specifications</u> Like the Master, the Delay Generator is housed in a two units ESONE plug-in (see fig. 1). The second Market Constants of A Output signal 1F' : - period equal to P2 of Master Generator - delay (with respect to 1F of Master Generator) : variable between 0 and the instantaneous value of P2, up to a maximum of 9999 /us), by steps of 1 /us (fig. 3). - instantaneous delay value displayed by 7 segment indicators. Output signal 6F' : - period equal to P1 of Master Generator, synchronous with 1F'. Common characteristics of both output signals : - pulse width : 2 μ s \pm 10 % March - quiescent level : + 15 V - pulses negative from + 15 V to ground

> capability of driving up to 10 MHTL circuits (I sink = 30 mA) from each output connector, including the test outputs
> slew rate : adjustable between 0.5 s and 1.5 s by front panel potentiometer

Clock frequency : - 1 Mhz, supplied by Master Generator

Operating modes :

- Manual : presetting and loading of the delay value by operator
- Automatic : the functions are computer controlled.

4) BIBLIOGRAPHY

 [1] Report EUR 4881e : A modular drive system for neutron choppers with digital speed and phase-control, by
 J. Eder and D. Röbbelen

[2] Report EUR to be published : Interface zur Steuerung einer Multichopperanlage by J. Eder and H. Brak.







- 19 -





.

- 21 -



- 22 -



- 23 -

'



- 24





- 26 -



NOTICE TO THE READER

All scientific and technical reports published by the Commission of the European Communities are announced in the monthly periodical "euro-abstracts". For subscription (1 year: BF.1025) or free specimen copies please write to:

> Office for Official Publications of the European Communities Boîte postale 1003 Luxembourg (Grand-Duchy of Luxembourg)

To disseminate knowledge is to disseminate prosperity — I mean general prosperity and not individual riches — and with prosperity disappears the greater part of the evil which is our heritage from darker times.

Alfred Nobel

SALES OFFICES

The Office for Official Publications sells all documents published by the Commission of the European Communities at the addresses listed below, at the price given on cover. When ordering, specify clearly the exact reference and the title of the document.

UNITED KINGDOM

H.M. Stationery Office P.O. Box 569 London S.E. 1 — Tel. 01-928 69 77, ext. 365

BELGIUM

Moniteur belge — Belgisch Staatsblad Rue de Louvain 40-42 — Leuvenseweg 40-42 1000 Bruxelles — 1000 Brussel — Tel. 12 00 26 CCP 50-80 — Postgiro 50-80

Agency :

Librairie européenne — Europese Boekhandel Rue de la Loi 244 — Wetstraat 244 1040 Bruxelles — 1040 Brussel

DENMARK

J.H. Schultz — Boghandel Møntergade 19 DK 1116 København K — Tel. 14 11 95

FRANCE

Service de vente en France des publications des Communautés européennes — Journal officiel 26, rue Desaix — 75 732 Paris - Cédex 15° Tel. (1) 306 51 00 — CCP Paris 23-96

GERMANY (FR)

Verlag Bundesanzeiger 5 Köln 1 — Postfach 108 006 Tel. (0221) 21 03 48 Telex: Anzeiger Bonn 08 882 595 Postscheckkonto 834 00 Köln

GRAND DUCHY OF LUXEMBOURG

Office for Official Publications of the European Communities Bolte postale 1003 — Luxembourg Tel. 4 79 41 — CCP 191-90 Compte courant bancaire: BIL 8-109/6003/200

IRELAND

Stationery Office — The Controller Beggar's Bush Dublin 4 — Tel. 6 54 01

ITALY

Libreria dello Stato Piazza G. Verdi 10 00198 Roma — Tel. (6) 85 08 CCP 1/2640

NETHERLANDS

Staatsdrukkerij- en uitgeverijbedrijf Christoffel Plantijnstraat 's-Gravenhage — Tel. (070) 81 45 11 Postgiro 42 53 00

UNITED STATES OF AMERICA

European Community Information Service 2100 M Street, N.W. Suite 707 Washington, D.C. 20 037 — Tel. 296 51 31

SWITZERLAND

Librairie Payot 6, rue Grenus 1211 Genève — Tel. 31 89 50 CCP 12-236 Genève

SWEDEN

Librairie C.E. Fritze 2. Fredsgatan Stockholm 16 Post Giro 193, Bank Giro 73/4015

SPAIN

Libreria Mundi-Prensa Castello 37 Madrid 1 — Tel. 275 51 31

OTHER COUNTRIES

122373981

Office for Official Publications of the European Communities Bolte postale 1003 — Luxembourg Tel. 4 79 41 — CCP 191-90 Compte courant bancaire: BIL 8-109/6003/300

OFFICE FOR OFFICIAL PUBLICATIONS OF THE EUROPEAN COMMUNITIES Boite postale 1003 — Luxembourg 6510

CDNA05112ENC